

PRO CLEAN-6



Working for Vinyl Records

Product Information

SURFACTANTS And Nothing Else

A better cleaning performance is now possible without alcohol. Special surfactants are the key.

Pro Clean-6 contains 3 of the latest generation of Lab Grade surfactants from the leading manufacturer of chemical reagents.

Wetting, detergent and emulsifier agents, all together.

Pro Clean-6 is a water based concentrate. It is alcohol free cleaning solution designed by Chemists and Record Collectors. It is suitable for all records including Vinyl, Shellac, Styrene, Lacquer and Acetates.

When mixed with distilled water the concentrate makes 6 litres of record cleaning fluid.

Features

- ✓ Removes a wide range of dirt, increasing bass response and improving soundstage clarity from LP playback.
- ✓ It does not leave any residues in the groove and the stylus thanks to its high purity formula without unnecessary compounds.
- ✓ Anti-static formula.
- ✓ Can clean hundreds of records very cost-effective.
- ✓ Fluid formula, clear, colourless and odourless. Mixed with distilled water it does not get cloudy and particles do not separate.

All Purpose Cleaning Methods

Pro Clean-6 is suitable for manual cleaning, bath cleaning, vacuum machines and Ultra-sonic cleaning machines.



Pro Clean-6 bottle

*Pro Clean-6 does not contain
brighteners, solvents, colourants,
fragrance or aggressive compounds*

Customer Feedback

"This stuff is magic, I have used a lot of cleaning fluids and this is the best"

"I'm impressed, it's far closed to alcohol than I expected for cleaning performance, works very well with styrene and vinyl records"

"I did try PC6 and I love it, might be the best I've used, very effective and pure"

Reviews

TNT-Audio
Review



The Essential Formula – Less is More

Pro Clean-6 contains the exact quantity of balanced compounds to allow residue-free cleaning.

The powerful surfactants just need an optimal quantity for a great job. This way the cleaner leave no trace on the record.

**No Alcohol
Safe and Effective**



Concentrate is made with biodegradable and non-toxic raw materials.

Water based formula increases Record longevity. Excellent alternative to conventional cleaners made with alcohol.

Product Specifications

<i>Dosage (for distilled water)</i>	<i>Ratio 1:100 for any amount. Examples: 10 ml. for 1 litre (1.000 ml.) of distilled water 1 ml. for 0,1 litres (100 ml.) of distilled water</i>
<i>Suitable for</i>	<i>Vinyl Records, including picture and coloured records Shellac recordings (78 rpm) Styrene, Lacquer and Acetate records</i>
<i>Cleaning methods</i>	<i>Vacuum machines (aspiration) Knosti type machines (bath machines) Ultra-Sonic cleaning machines Manual cleaning, including immersion</i>
<i>Can clean</i>	<i>Vacuum Machines: 3,0 ml. per side = 1.000 records Manual washing (spray): 1 ml. per side = 3.000 records Knosti and Ultrasounds depends on the usage</i>

Sales Specifications

<i>Bottle</i>	<i>60 ml. PEHD square bottle Polypropylene cap</i>
<i>Box</i>	<i>1 Syringe of 5 ml. 1 Syringe of 2 ml. 1 Dispensing needles Instructions Guide (EN, FR, DE, ES)</i>
<i>Measures</i>	<i>6 x 6 x 9 cm (L x W x H)</i>
<i>Gross weight</i>	<i>110 grams</i>
<i>SKU</i>	<i>WPC-6</i>
<i>EAN</i>	<i>0763846976365</i>

Winyl Record Cleaner
winylcleaner.com

Contact

contact@winylcleaner.com
winylcleaner.com